

Specification

No.	Item	Description	UOM	Quantity	Picture (for reference)
1	a. Anatomical model for Implant	<p>MODELARM: Subcutaneous Implant Training Arm</p> <p>Reproductive Implant Training Arm A compact simulator designed to insert and remove Levonorgestrel (Norplant) implants.</p> <p>Reproductive Implant Training Arm Features:</p> <ul style="list-style-type: none"> Upper left arm on base Soft foam insert simulates soft arm tissue Soft foam insert can be rotated 360° allowing multiple insertion exercises 5 tubular inserts One extra latex skin Instructions 	Set	8	
	b. Extra Latex Skin (refers to item no. 1.a)	Specification is the same as extra latex skin under the package/set of item no 1.a	Pcs	20	
2	Anatomical model for IUD	<p>Desk Top I.U.D. Trainer - Female Pelvic Organs I</p> <p>Desk Top I.U.D. Trainer - Female Pelvic Organs I - Medium Skin Tone This compact Desktop Trainer Female Pelvic Organs Model illustrates the position the IUD assumes once inserted correctly into the uterus.</p> <p>Desk Top Trainer Female Pelvic Organs I Features:</p> <ul style="list-style-type: none"> Coronal section of uterus ovaries and fimbriae Uterus covered by clear plastic window to permit easy viewing of correct IUD placement Lightweight and portable Instruction manual Included 	Set	8	
3	Anatomical model for IUD/Minilaparotomy	<p>Gynecologic Simulator</p> <p>The anatomical model should has ability to demonstrate multiple gynecologic procedures as well as practice laparoscopic examination and minilaparotomy. The simulator should be Latex-free.</p> <p>The complete package should include all uterus models and cervix.</p>	Set	8	
4	a. Anatomical model for Vasectomi	<p>No Scalpel Vasectomy (NSV) Simulator</p> <p>NSV - Non-Scalpel Vasectomy Simulator Model Features:</p> <ul style="list-style-type: none"> Two removable scrotal skins Two testicles Two long vas assemblies (vas deferens) can be advanced as needed for NSV exercises 	Set	8	
	b. Extra scrotal skins (refers to item no 4.a)	Specification is the same as extra scrotal skins under the package/set of item no 4.a	Pcs	20	
	c. Extra vas deferens/long vas assembly (refers to item no 4.a)	Specification is the same as vas deferens/long vas assemblies under the package/set of item no 4.a	Pcs	20	
5	Anatomical model for IUD	<p>IUD Insertion Trainer Uterus</p> <p>IUD Insertion Trainer Uterus Features:</p> <ul style="list-style-type: none"> Normal uterus Clear plastic window permits easy viewing of IUD Plastic window tilts open to permit removal of IUD 	Set	8	